

for the COVVI Hand - Friction Wrist

We here with declare, under the manufacturer's sole responsibility, that the products described in this document have been evaluated against **The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012 (Statutory Instrument 2012 No. 3032)**. These regulations implement the restriction of the use of certain hazardous substances in electrical and electronic equipment placed on the market in Great Britain. The compliance status of the product is confirmed by the sections below.

PRODUCT NAME:	COVVI Hand - Friction Wrist
LEGAL MANUFACTURER:	COVVI LTD, Unit 4 (Direct House), Quayside Business Park, George Mann Road, Leeds, LS10 1DJ, United Kingdom
CONFIGURATIONS & VARIANTS:	As per Appendix I - Product Listing/Schedule
INTENDED PURPOSE:	The COVVI Hand - Friction Wrist is an external, motor-driven prosthetic hand with a Friction Wrist designed to assist individuals with unilateral or bilateral upper limb deficiencies. Specifically crafted with anatomically correct dimensions, the COVVI Hand - Friction Wrist facilitates activities of daily living through controlled movement of programmed grip patterns.

Based on the information provided by our supply lines and our certain knowledge pertaining to our own processes, all products described in this document conform to the requirements of The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012 (Statutory Instrument 2012 No. 3032) and do not contain any of the below-mentioned 10 Hazardous substances above the specified limits.

Restricted substances	Maximum Threshold
Cadmium	100 ppm (0.01 weight%)
Mercury	1000 ppm (0.1 weight %)
Hexavalent chromium	1000 ppm (0.1 weight %)
Lead	1000 ppm (0.1 weight %)
Polybrominated biphenyls (PBB)	1000 ppm (0.1 weight %)
Polybrominated diphenyl ethers (PBDE)	1000 ppm (0.1 weight %)
Bis(2-ethylhexyl phthalate (DEHP)	1000 ppm (0.1 weight %)
Benzyl butyl phthalate	1000 ppm (0.1 weight %)
Dibutyl phthalate (DBP)	1000 ppm (0.1 weight %)
Di isobutyl phthalate (DIBP)	1000 ppm (0.1 weight %)

Name Simon Pollard **Position** Group CEO



Signed _____ **Date** 12-Dec-2025 **Place** Leeds, United Kingdom

Who is the natural and legal person with responsibility for the design, manufacture, packaging and labelling before the device is placed on the market under this manufacturer's name regardless of whether these operations are carried out by the manufacturer or on his behalf by a third party.

Our statements in this document regarding RoHS compliance do not extend to, or apply to any product subjected to unintended contamination, misuse, neglect, or accident.

COVVI

COVVI House, 4 Quayside Business Park, George Mann Road, Leeds, LS10 1DJ
Switchboard: +44 (0)203 9499 500 Website: www.covvi.com Website: www.covvi-robotics.com

Appendix I - Product Listing/Schedule (Configurations and Variants)

Product Code	Description
CVSLSFBLRG	Friction Wrist Small Left Hand - Black Glove, Rose Gold Cover
CVSLSFBLCA	Friction Wrist Small Left Hand - Black Glove, Carbon Cover
CVSLSFBLTG	Friction Wrist Small Left Hand - Black Glove, Titan Grey Cover
CVSLSFBLWH	Friction Wrist Small Left Hand - Black Glove, White Cover
CVSRSFBLRG	Friction Wrist Small Right Hand - Black Glove, Rose Gold Cover
CVSRSFBLCA	Friction Wrist Small Right Hand - Black Glove, Carbon Cover
CVSRSFBLTG	Friction Wrist Small Right Hand - Black Glove, Titan Grey Cover
CVSRSFBLWH	Friction Wrist Small Right Hand - Black Glove, White Cover
CVSLMFBLRG	Friction Wrist Medium Left Hand - Black Glove, Rose Gold Cover
CVSLMFBLCA	Friction Wrist Medium Left Hand - Black Glove, Carbon Cover
CVSLMFBLTG	Friction Wrist Medium Left Hand - Black Glove, Titan Grey Cover
CVSLMFBLWH	Friction Wrist Medium Left Hand - Black Glove, White Cover
CVSRMFBLRG	Friction Wrist Medium Right Hand - Black Glove, Rose Gold Cover
CVSRMFBLCA	Friction Wrist Medium Right Hand - Black Glove, Carbon Cover
CVSRMFBLTG	Friction Wrist Medium Right Hand - Black Glove, Titan Grey Cover
CVSRMFBLWH	Friction Wrist Medium Right Hand - Black Glove, White Cover
CVSLLFBLCA	Friction Wrist Large Left Hand - Black Glove, Carbon Cover
CVSLLFBLTG	Friction Wrist Large Left Hand - Black Glove, Titan Grey Cover
CVSLLFBLWH	Friction Wrist Large Left Hand - Black Glove, White Cover
CVSRLFBLCA	Friction Wrist Large Right Hand - Black Glove, Carbon Cover
CVSRLFBLTG	Friction Wrist Large Right Hand - Black Glove, Titan Grey Cover
CVSRLFBLWH	Friction Wrist Large Right Hand - Black Glove, White Cover
LNSLSF	Loan Friction Wrist Left Small Hand
LNSRSF	Loan Friction Wrist Right Small Hand
LNSLMF	Loan Friction Wrist Left Medium Hand
LNSRMF	Loan Friction Wrist Right Medium Hand
LNSLLF	Loan Friction Wrist Left Large Hand
LNSRLF	Loan Friction Wrist Right Large Hand
Accessories	
CFWLR	COVVI Friction Wrist Lamination Ring
CFWSB	COVVI Friction Wrist Switch Block
LFWSB	Loan Friction Wrist Switch Block
CFWIT	COVVI Friction Wrist Installation Tool